



**Center for Clinical and
Translational Research
Standard Operating Procedure**



Section: **Clinical Research Center**

Date Created: **August 2nd, 2021**

Title: **Eyewash Station Check**

Version Date: **January 1, 2023**

SOP Number: **CO22**

PURPOSE: The purpose of this standard operating procedure (SOP) is to describe the process for weekly eyewash station inspections to verify operation and ensure that water is available and in the correct temperature range.

SCOPE: This SOP applies to all staff that are inspecting the eyewash station in the CRC.

PERSONNEL RESPONSIBLE: Research Associates

PROCEDURES:

1. Ensure access to eyewash is unobstructed.
2. Verify that eyewash station is clean and ready for use.
3. Turn on eyewash station and allow water to flow for 15 minutes.
4. Take the temperature of the water after 15 minutes. The water temperature should be between 60° and 100°F.
5. Document the inspection by putting the date, temperature, and your initials. This inspection should be recorded weekly.
6. Report problems to Facilities at 402-552-3347.

RESOURCES:

ANSI standard Z358.1-2014

Department Approval

Signed *Serena Gaires*
Research Nurse Manager

Signed: 2/9/2023

Signed *[Signature]*
Assistant Vice Chancellor for Clinical Research

Signed: 2/7/2023