

PREPARATION FOR ELUTRIATION (April 2020)

1. Assemble chambers and 3 transfer tubes. Follow diagram, starting at the bottom. Lubricate O-rings using Beckman vacuum grease silicon.
2. Ready the centrifuge. Turn on. The settings should remain from previous time. Time = hold (H); Temp = 18° – 25 °C; RPM = 1960. Attach chambers to rotor and secure cable to rotor. Connecting tubing: Out = top; In: side.
3. Clean system tubing: Replace stopcock and syringe. Place new 1 ml aspirating pipette in a bottle of 95% EtOH. Make sure pipette is all the way down to the bottom of bottle but not touching it. Replace new short 1 ml pipette at the outlet. Turn on pump (forward) and flush with EtOH with speed setting @ 2. You should see flow into waste collection. Invert chamber, tap to free any bubbles.
4. After ~200 ml of EtOH has run through the system, turn off the pump. Replace EtOH with 1L bottle of PBS (always keep this bottle filled). Turn on the pump (set @ 2). Pour PBS into the syringe and elute this to rid any air bubbles that might be in the stopcock. Close the stopcock before all of the PBS has drained and continue to run PBS through the system. Make sure you have enough PBS on hand!
5. Continue the PBS flush for ~400 ml (speed setting @2). During this flush, run centrifuge at 1960 rpm, 18 °C (recall program #1), pinching tubing continuously to eliminate any trapped air bubbles in the tubing. When the speed gets up to 1960 rpm, stop the centrifuge, and turn the pump up to 4.0 to make sure it can handle the pressure. Turn the pump back down to 1.0.
6. Check chamber for any trapped air bubble. Repeat step 5 to rid of air bubbles. Turn off the pump until ready for elutriation. (Pump should always be on before and after centrifuge operation)
7. Turn on Coulter counter. Empty waste if more than half full and fill Isoton chamber 2/3 full.
8. Flush electrode chamber with Isoton in cup (10 ml). Press START to count. Count should be ~500.
9. To count cells, add 10 ml Isoton to 20 µl cells. Place under electrode. Press START to count. Press OUTPUT to read count number.

$$\text{Cell concentration (cells/ml)} = \text{raw count} \times 10^3 \times 10 \text{ (dilution factor)}$$

To Count Cells

Add 10 ml Isoton to 20 μ l cells. Place under electrode. Press **START** to count. Press **OUTPUT** to see cell count and profile.

CLEANING AFTER ELUTRIATION

CENTRIFUGE

1. Run 400 ml PBS through the system with pump set at 2.
2. Run 200 ml EtOH afterward.
3. Disassemble the chamber, clean the chamber by first soak in 10% Solution 555 for 10 min. then rinse off with water and EtOH at the end.
4. Clean all others with EtOH
5. Turn off centrifuge
6. Clean hood with EtOH.
7. **Empty vacuum flask**, make sure it was bleached before going down the sink. Rinse and fill with bleach and assemble.

Coulter Counter

1. Count 2x with Coulter Clenz. Leave cup in there.
2. Turn all three instruments off.
3. Wipe surface off with EtOH.